



Piccolo is perfect for everyone, removing plaque and debris for a thorough clean and supports gum health.

Piccolo is also, especially good for patients with limited dexterity and those with dental prosthetics

---

## Join our user evaluation



User Evaluation  
Landing page



PICCOLO

## Follow us



PICCOLO



PICCOLO

## Contact Us

Trevor Martin  
Email: [tmartin@tt3.com.au](mailto:tmartin@tt3.com.au)  
Ph: 0407512378  
Address: TT3  
C/O The Nest,  
Level 2, 180 King Street ,  
Newcastle, 2300,  
NSW, Australia

Meeting Request



PICCOLO

# PICCOLO



A world first in  
interdental cleaning

## Powered

5 Easy steps

1. Turn the Piccolo on
2. Position the tip at the interdental space on the gum line
3. Squeeze the grip
4. Release the grip
5. Move to the next interdental space. No mess, no pain, quiet, fast



## Why Piccolo

### Perfect for everyone

Including those with, tight spaces between their teeth, orthodontic appliances, implants, crowns and periodontal conditions

### Universal tip fits all

### Thorough

Removes plaque & food debris effortlessly and supports gum health




### Gentle & Safe

Piccolos' tip with smart motion means no more, black triangles from over brushing or using too much force, or micro abrasion of tooth surfaces compared to wood picks, oral irrigators, floss or interdental brushes. And the user has full control over the depth of insertion

## Tip Design

Piccolo's tapered helical tip, made from, food grade recyclable, thermoplastic elastomer works with **SMART MOTION** technology, combining 1500rpm rotational and 12 mm axial movement' to clean thoroughly and gently while enhancing user comfort.



### Affordable and ecofriendly

Rechargeable battery and only one tip means lower cost than leading oral irrigators and over a year less expensive than interdental brushes. Less waste

### Convenient

One-handed use, no mess, no stress, fast and easy and easy to travel with